



TECHNICAL SPECIFICATIONS

Seat Pad Size: 17.5" / 445 mm wide x 15" / 381 mm deep

Armrests: adjustable to 21.75" / 552 mm between arms

Footrest Size: 10" / 254 mm long x 9" / 229 mm wide



No matter what the weather, Bruno's Elite Outdoor Curved stairlift helps you enjoy the freedom of the great outdoors. Custom-crafted to precisely follow the angles of your stairs, the Elite offers Bruno's renowned reliability and smooth, safe ride. The weather-protected design ensures top-notch performance year round.

- Weight capacity: 400 lb/ 181 kg
- Performance tested for 0 F - 125 F (-18 C to 52 C)
- Durable outdoor paint protects all steel surfaces
- Marine-grade vinyl padded seat and armrests
- Weather-resistant cover protects the unit from the elements when stairlift not in use
- Flip up arms, seat, and footrest create plenty of extra space on steps
- Obstruction safety sensors, seat belt ensure safety
- Two 12v batteries ensure dependable performance even in a power outage
- Soft start/stop gives smooth ride from start to finish
- Plush, generous size seat with multiple height adjustments
- Onboard audio/visual diagnostics for easy service
- Offset swivel seat makes getting on/off safe and easy
- Gold warranty: five year major components; two year parts

See reverse side for additional features. Learn more at www.bruno.com

ELITE OUTDOOR CURVED FEATURES



Footrest/carriage safety sensors stop the unit when obstruction encountered.



Offset swivel seat makes getting on/off easy.



Flip up arms, seat and footrest create plenty of extra space on steps.



Weather-resistant cover provides durable protection.



Two wireless call/send controls.



Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066
www.bruno.com

© Bruno Independent Living Aids, Inc. 2015

MKT-SL-5