



TECHNICAL SPECIFICATIONS

Standard Seat Pad Size: 17.5" / 445 mm wide x 15" / 381 mm deep

Larger Seat Pad Option: 21.5" / 546 mm wide x 16" / 406 mm deep

Standard Footrest Size: 10" / 254 mm long x 9" / 229 mm wide

Larger Footrest Option: 12" / 305 mm long x 11" / 279 mm wide



When you want the highest quality stairlift combined with a stylized appearance, you want the Bruno Elite. The Elite offers worry-free access to all levels of your home. And the Elite's stylized design blends with your home's environment. Consider custom upholstery the chair in your own fabric.

- Weight capacity: 400 lb/ 181 kg
- Plush, generous-size seat with multiple height adjustments
- Adjustable width armrests to 21.75" / 552 mm between arms
- Offset swivel seat makes getting on/off safe and easy
- Under-armrest rocker switch for easy operation
- Anodized rail with covered gear rack gives a sleek appearance
- Flip up arms, seat, and footrest create plenty of extra space on steps
- Adjustable footrest for maximum comfort
- Tight rail installation: 2" / 51 mm wide rail installs to within 4" / 101 mm of the wall
- Continuous charge strips along rail power unit wherever it stops
- Two wireless call/sends allow unit to be controlled remotely
- Two 12v batteries ensure dependable performance even in a power outage
- Soft start/stop gives a smooth ride from start to finish
- Onboard audio/visual diagnostics for easy service
- Obstruction safety sensors, seat belt ensure safety
- Gold warranty: five year major components; two year parts

See reverse side for optional features. Learn more at www.bruno.com


ELITE STRAIGHT



ELITE STRAIGHT OPTIONS

Standard tan vinyl with five optional upholstery selections. Or, supply your own fabric and Bruno will custom upholster the chair for you.

- Wider seat
- Wider footrest



Power folding footrest automatically flips up/down when seat is raised/lowered.



Power or manual folding rails for narrow hallway or when doorway is at the bottom of stairs. Manual operation or push button automatic.



Power swivel seat for effortless exit. Controlled on chair arm or wireless call/send.



Powder-coated champagne painted rail.



Bruno Independent Living Aids, Inc.
 1780 Executive Drive, Oconomowoc, WI 53066
www.bruno.com

© Bruno Independent Living Aids, Inc. 2015

MKT-SL-2